

# The AutoBio Instruction Set Manual

2022

# Preface

This is version 0.1 (draft) of the document describing the AutoBio Instruction Set Architecture (AutoBio ISA). The purpose of the Instruction Set is to provide an interface between humans and machines in the field of biology, i.e., to formulate a clearly defined exchangeable format for biology experiments.

Our aim is to use AutoBio ISA as an alternative to biology protocols in natural languages, and this instruction set is specially designed for experiment standardization.

Currently, this instruction set is far from finished, refinement is welcomed. For more information, please visit

<https://www.overleaf.com/6629364199jbrwrrxyyrg>

# Inhaltsverzeichnis

<b>1</b>	<b>Instruction Overview</b>	<b>5</b>
1.1	AutoBio ISA Overview . . . . .	5
1.2	Instruction Encoding . . . . .	6
1.3	Reagent Encoding . . . . .	7
<b>2</b>	<b>Data Types</b>	<b>8</b>
2.1	Basic Types . . . . .	8
2.1.1	Booleans . . . . .	8
2.1.2	Numbers . . . . .	8
2.1.3	Temperature Specifications . . . . .	8
2.1.4	Pressure Specifications . . . . .	8
2.1.5	Volume Specifications . . . . .	8
2.1.6	Duration Specifications . . . . .	8
2.2	Structural Types . . . . .	8
2.2.1	Reagents . . . . .	8
2.2.2	Conditions . . . . .	8
<b>3</b>	<b>Basic Actions</b>	<b>9</b>
3.1	add . . . . .	10
3.2	concentrate . . . . .	11
3.3	incubate . . . . .	12
3.4	elute . . . . .	13
3.5	wash . . . . .	14
3.6	transfer . . . . .	15
3.7	spin . . . . .	16

*Inhaltsverzeichnis*

3.8	collect . . . . .	17
3.9	shake . . . . .	18
3.10	suspend . . . . .	19
3.11	dialyze . . . . .	20
<b>4</b>	<b>Pseudo Actions</b>	<b>21</b>
4.1	repeat . . . . .	21
<b>5</b>	<b>Analytical Actions</b>	<b>23</b>

# 1 Instruction Overview

The AutoBio-ISA is a new instruction set architecture (ISA) designed to support computer-aided biological experimental automation. Our goals in defining AutoBio-ISA include the following:

- A *real* ISA suitable for biological research, not just limited to computer science.
- An ISA that avoids “*over-natural*” in the current protocol standard. We aim to make each step as clear as possible.
- An ISA that *omits* the execution details, e.g., the specified well for washing.
- An ISA which exploits the parallel nature in the art of experiment.

## 1.1 AutoBio ISA Overview

The central part of the base AutoBio-ISA is the *execution graph* and *instruction*, we split the entire protocol into different semantic blocks, and translate the procedures written in natural language inside the block into AutoBio instructions. We exploit the block-level parallelism by checking the mutual interference between each block, creating the dependency graph, and finding the critical path on this graph, see chapter ??.

## 1 Instruction Overview




Abbildung 1.1: An example of instruction flow

## 1.2 Instruction Encoding

An *instruction* describes a specific form of a laboratory operation, including the operation itself (called a “*mnemonic*”), the input of the operation (called the “*slot*” of the action), and the output of the operation (called the “*emit*” of the action). As an example, the following instruction describes the procedure of reagent addition followed by incubation. The original form of the instruction flow is represented in a graph form where the edges denote the reagent transferring from one step to another.

As shown in Abbildung 1.1, the only constraint between instructions is the flow of reagent. An instruction emits reagents, and another instruction catches them in its slots, forming a flow edge. To serialize the graph into a sequence, we assume that each instruction has exactly *one* or no output, and assign a unique number to each instruction. The sequential version of the instruction is

```
15: ADD %14, %5
16: INCUBATE %15
17: PLACE %16 %4 "freeze"
```

where the instruction ADD has the number 15 and its input slot catches the reagents emitted from the instruction 14 and 5. After incubation, the resulting reagent is placed in the condition

## 1 *Instruction Overview*

specified in 4. The optional field "freeze" comes from the protocol text, indicating the PLACE operation in the instruction list is originally described as "freeze".

The instruction has the form MNEMONIC SLOTS LITERAL, where the optional field LITERAL records the verb used in the original protocol, the MNEMONIC field is a pre-defined value selected from a set of mnemonic we defined in Chapter 3.

### 1.3 Reagent Encoding

TODO

## **2 Data Types**

This chapter describes the data types defined for AutoBio-ISA. There are basic types and structural types.

### **2.1 Basic Types**

#### **2.1.1 Booleans**

#### **2.1.2 Numbers**

#### **2.1.3 Temperature Specifications**

#### **2.1.4 Pressure Specifications**

#### **2.1.5 Volume Specifications**

#### **2.1.6 Duration Specifications**

### **2.2 Structural Types**

#### **2.2.1 Reagents**

#### **2.2.2 Conditions**



## 3 Basic Actions

This chapter describes the basic instructions used in biological experiments. We aim to reduce concrete actions into abstract operations, without loss of their executability.

An instruction may have some slots (inputs) or emits (outputs), and each slot or emit has a specified data type attached to it.

The REQ flag denotes that this argument is *Required*, the OPT flag denote that this argument is *Optional*, and the RET flag denote that this argument is an emit (return value).

The REG, COND, BOOL flag represents the type of *Reagent*, *Condition* and *Boolean*, respectively.

Each instruction corresponds to an instruction node in the execution graph.

## 3.1 add

```
r3 = ADD r1, r2, Ordered=True
```

### Description

Add a reagent to another reagent. This action is the basic abstraction of adding a reagent to another. Examples include mixing reagents, applying a drop, etc.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The <i>target</i> reagent.
r2	Slot	REQ REG	The <i>additive</i> reagent.
b3	Slot	OPT BOOL	r1 and r2 are interchangeable.
r4	Emit	RET REG	The <i>result</i> reagent.

### Examples

P/I	Content
Proto.	Apply a drop of normal goat serum to each specimen.
Instr.	ADD [specimen], [goat serum], b3=True <sup>1</sup>
Proto.	Load supernatant on to 3 ml equilibrated Ni-NTA column.
Instr.	ADD [Ni-NTA], [supernatant], b3=True, LOAD

<sup>1</sup>The [specimen] and [goat serum] are two reagent objects holding their solutes, volumes, etc..

## 3.2 concentrate

r3 = CONCENTRATE r1, c2

### Description

Concentrate a reagent to a specified condition.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be concentrated.
r2	Slot	REQ COND	The target condition.
r3	Emit	RET REG	The <i>result</i> reagent.

### Examples

P/I	Content
Proto.	Adjust the cells concentration into $1.5 \times 10^5$ cells/mL.
Instr.	CONCENTRATE [cells], [concentration= $1.5 \times 10^5$ ]
Proto.	Concentrate eluate in speedvac to 22.5 $\mu$ L
Instr.	CONCENTRATE [eluate], [volume=22.5 $\mu$ L]

### 3.3 incubate

r3 = INCUBATE r1, c2

#### Description

Incubate a reagent in a specified environment.

#### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be incubated.
r2	Slot	REQ COND	The incubation environment.
r3	Emit	RET REG	The <i>result</i> reagent.

#### Examples

P/I	Content
Proto.	Incubate at 37°C for 2-3 minutes.
Instr.	INCUBATE [cells], [temp=37C, dur=2-3min]
Proto.	Incubate the mixture for 5 min at room temperature.
Instr.	INCUBATE [cells], [temp=ROOM, dur=5min]

## 3.4 elute

`r4 = ELUTE r1, r2, c3`

### Description

Elute the reagent `r1` in the eluent `r2` under the condition `c3`, the resulting reagent is given in `r4`.

### Slots & Emits

Argument	E/S	Type	Description
<code>r1</code>	Slot	REQ REG	The reagent to be eluted.
<code>r2</code>	Slot	REQ REG	The eluent reagent.
<code>c3</code>	Slot	OPT COND	The environment condition for elution.
<code>r3</code>	Emit	RET REG	The result eluate reagent.

### Examples

P/I	Content
Proto.	Elute the histone complex with 250 $\mu$ L freshly prepared elution buffer.
Instr.	ELUTE [complex], [elution buffer]
Proto.	Elute RNA in 20 $\mu$ L of Elution Solution A.
Instr.	ELUTE [rna], [elution solution A]

## 3.5 wash

`r4 = WASH r1, r2, c3`

### Description

Wash the reagent `r1` in the fluid reagent `r2` under the condition `c3`, the resulting reagent is given in `r4`.

### Slots & Emits

Argument	E/S	Type	Description
<code>r1</code>	Slot	REQ REG	The reagent to be washed.
<code>r2</code>	Slot	REQ REG	The wash fluid.
<code>c3</code>	Slot	OPT COND	The environment condition for washing.
<code>r3</code>	Emit	RET REG	The washed reagent.

### Examples

P/I	Content
Proto.	Wash the cell in 1ml PCR-grade water.
Instr.	<code>WASH [cell], [PCR-grade water]</code>
Proto.	Wash the digested tissues with HBSS containing 5 mM EDTA.
Instr.	<code>WASH [tissues], [HBSS]</code>

## 3.6 transfer

```
r3 = TRANSFER r1, r2
```

### Description

Transfer the reagent r1 to the reagent r2, the resulting reagent is given in r3.<sup>2</sup>

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be transferred.
r2	Slot	REQ REG	The target container.
r3	Emit	RET REG	The transferred reagent.

### Examples

P/I	Content
Proto.	Transfer the aqueous phase to a tube containing 500 $\mu$ L isopropanol.
Instr.	TRANSFER [aqueous phase], [tube]
Proto.	Transfer the mixture to a phase-lock tube.
Instr.	TRANSFER [mixture], [phase-lock tube]

<sup>2</sup>The reagent r2 might be an empty container.

## 3.7 spin

r3 = SPIN r1, c2

### Description

Spin (or centrifugate) the reagent r1 under the condition c2, the resulting reagent is given in r3.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be centrifugated.
c2	Slot	REQ COND	The centrifugation settings.
r3	Emit	RET REG	The resulting reagent.

### Examples

P/I	Content
Proto.	Centrifuge for 5 min 30 sec at 230 rcf to pellet cells.
Instr.	SPIN [cells], [5m30s,230rcf]
Proto.	Spin down the plates at 110 rcf for 6 min at room temperature.
Instr.	SPIN [plate], [6m,110rcf,room temp]

### Important Notes

We treat *Spin* and *Centrifuge* as the same operation.



## 3.8 collect

```
r3 = COLLECT r1, r2
```

### Description

Collect the reagent `r1` to the reagent (or container) `r2`, the resulting reagent is given in `r3`.

### Slots & Emits

Argument	E/S	Type	Description
<code>r1</code>	Slot	REQ REG	The reagent to be collected ('collect-from').
<code>r2</code>	Slot	OPT REG	The target reagent or container ('collect-to').
<code>r3</code>	Emit	RET REG	The resulting reagent.

### Examples

P/I	Content
Proto.	Collect suspensions into 2 x 250 mL conical bottles.
Instr.	COLLECT [suspension], [conical bottle]
Proto.	Collect resulting emulsion via PE2 tubing into a 50 mL Falcon tube.
Instr.	COLLECT [emulsion], [Falcon tube]
Proto.	Collect livers from e13.5 mouse embryos.
Instr.	COLLECT [e13.5 mouse embryos], [livers]

### Important Notes

Clarification between `COLLECT` and `TRANSFER` needs to be provided. There exists ambiguity in the above examples.

## 3.9 shake

r3 = SHAKE r1, c2

### Description

Shake the reagent r1 under the condition c2, the resulting reagent is given in r3.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be shaken.
c2	Slot	REQ COND	The shake settings.
r3	Emit	RET REG	The resulting reagent.

### Examples

P/I	Content
Proto.	Shake the lymph ducts gently.
Instr.	SHAKE [lymph ducts], [gently]
Proto.	Shake vigorously for 15 seconds.
Instr.	SHAKE [cell], [vigorously,15sec]

### Important Notes

SHAKE might be both an action or a field of COND (with shaking).

## 3.10 suspend

r3 = SUSPEND r1, r2

### Description

Suspend the reagent r1 using the reagent r2, the resulting reagent is given in r3.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be suspended.
r2	Slot	REQ REG	The suspension liquid.
r3	Emit	RET REG	The resulting reagent.

### Examples

P/I	Content
Proto.	Aspirate supernatant and suspend cells in PGC medium.
Instr.	SUSPEND [cells], [PGC medium]
Proto.	Suspend frozen yeast cell pellets in 5 ml lysis solution.
Instr.	SUSPEND [cells], [lysis solution]

## 3.11 dialyze

r3 = DIALYZE r1, r2, c3

### Description

Dialyze the reagent r1 using the reagent r2, the condition for dialysis is given in c3.

### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be dialyzed.
r2	Slot	REQ REG	The dialysis solution.
c3	Slot	OPT COND	The condition for dialysis.
r4	Emit	RET REG	The result reagent.

### Examples

P/I	Content
Proto.	Dialyze flow through in FPLC Buffer A.
Instr.	DIALYZE [flow], [FPLC Buffer A]
Proto.	Dialyze protein into appropriate experimental buffer.
Instr.	DIALYZE [protain], [appropriate buffer]

## 4 Pseudo Actions

### 4.1 repeat

r3 = REPEAT s1, r2

#### Description

Repeat the

#### Slots & Emits

Argument	E/S	Type	Description
r1	Slot	REQ REG	The reagent to be collected ('collect-from').
r2	Slot	OPT REG	The target reagent or container ('collect-to').
r3	Emit	RET REG	The resulting reagent.

#### Examples

#### Important Notes

Clarification between COLLECT and TRANSFER needs to be provided.

There exists ambiguity in the above examples.

#### 4 Pseudo Actions

P/I	Content
Proto.	Collect suspensions into 2 x 250 mL conical bottles.
Instr.	COLLECT [suspension], [conical bottle]
Proto.	Collect resulting emulsion via PE2 tubing into a 50 mL Falcon tube.
Instr.	COLLECT [emulsion], [Falcon tube]
Proto.	Collect livers from e13.5 mouse embryos.
Instr.	COLLECT [e13.5 mouse embryos], [livers]

## 5 Analytical Actions

measure

analyse

plot

fit